CanGaroo®

ENVELOPE INSTRUCTIONS FOR USE

Manufacturer:
Aziyo Biologics, Inc.
1100 Old Ellis Road, Suite 1200
Roswell, Georgia 30076 USA
Phone: +1 470-514-4080
Fax: +1 678-680-5486
www.aziyo.com

Authorized Representative in the
European Community:
ICON (LR) Limited
South County Business Park
Leopardstown, Dublin 18, D18 X5R3
Ireland

CE 0344

INDICATIONS FOR USE

The CanGaroo® Envelope is intended to securely hold a cardiac implantable electronic device to create a stable environment when implanted in the body. The devices that may be used with the CanGaroo® Envelope include pacemaker pulse generators, defibrillators, or other cardiac implantable electronic devices.

CONTENTS:

One (1) sterile, non-pyrogenic CanGaroo® Envelope

PRODUCT DESCRIPTION

The CanGaroo® Envelope is constructed from perforated, multilaminar sheets (4-ply) of decellularized, non-crosslinked, lyophilized extracellular matrix (ECM) derived from porcine small intestinal submucosa. The perforations in the ECM material allow the exit of any exudate. The ECM is assembled into pouch form with 5-0 polydioxanone (PDS) suture.

Tissue ingrowth into a properly sized CanGaroo® Envelope will promote stabilization of the device when compared to implantation with only standard fixation methods (e.g., sutures through the CIED header) or implantation with no fixation methods. By containing the CIED in the CanGaroo® and through tissue ingrowth, a stable environment is created in which there is less room for movement of the CIED in the tissue pocket. By creating a stable environment, the likelihood for complications associated with CIED movement are reduced. Such complications may include migration of the CIED, erosion of the CIED through the skin, or complications associated with Twiddler's syndrome.

HOW SUPPLIED

The CanGaroo® Envelope is supplied **STERILE**. Provided that the integrity of the sterile pouch is not compromised in any way, it serves as an effective sterile barrier until the "Use By" (expiration) date printed on the pouch.

The CanGaroo® Envelope is provided in the sizes listed in Table 1.

Table 1 – CanGaroo® Envelope Sizes

Size	Envelope Height (cm)	Envelope Width (cm)
SMALL (S)	5.0	5.4
MEDIUM (M)	6.5	6.9
LARGE (L)	8.0	6.9
EXTRA LARGE (XL)	9.5	6.9
EXTRA EXTRA LARGE (XXL)	8.9	10.8

PRECLINICAL DATA

In the rabbit model, the CanGaroo® Envelope demonstrated effectiveness in providing a barrier surrounding a CIED compared to a pacemaker canister alone in rabbits. Ingrowth of vascularized tissue was observed in the CanGaroo® Envelopes after six weeks. The ingrowth of tissue provided additional stabilization for the CIED canister in the rabbit model compared to CIEDs implanted without the CanGaroo® Envelope.

CONTRAINDICATIONS

- The CanGaroo® Envelope is derived from a porcine source and should not be used in patients with a known sensitivity to porcine material.
- 20ml gentamicin 40mg/ml should not be used as a hydration fluid for the CanGaroo® Envelope in patients with known hypersensitivity to gentamicin. A history of hypersensitivity or serious toxic reactions to other aminoglycosides may contraindicate use of gentamicin because of known cross-sensitivity of patients to drugs in this class.

WARNINGS AND PRECAUTIONS

- Only physicians qualified in the placement of electronic implantable devices should use this device.
- For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization will compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death.
- Discard all open or unused product.
- Device is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Do not use if storage conditions have not been maintained.
- The device must be used prior to the expiration date.
- Discard device if mishandling has caused possible damage or contamination.
- Always handle the device using aseptic technique.
- Do not use the CanGaroo® Envelope with fibrin glue or glue containing glutaraldehyde or other cross-linking agents.
- The use of the CanGaroo® Envelope has not been studied in the pediatric population. As with adult patients, an appropriately-sized CanGaroo® Envelope should be selected according to the size of the implantable device.
- Ensure that the device is properly hydrated prior to suturing. Without proper hydration, the ECM will tear and may not properly retain sutures.
- Once hydrated, the CanGaroo® Envelope should either be used or discarded. The CanGaroo® Envelope should not be rehydrated and reused.
- If 20ml gentamicin 40mg/ml is selected for hydration of the CanGaroo® Envelope, the contraindications, warnings, and precautions applicable to the use of gentamicin should be followed.
- If 20ml gentamicin 40mg/ml is selected for hydration of the CanGaroo® Envelope, it must be hydrated using aseptic technique by completely immersing it in the bowl of gentamicin for **2 minutes** prior to use. Hydrating the Envelope for less than 2 minutes may result in incomplete hydration and the Envelope may tear; hydration longer than 2 minutes may make the Envelope more difficult to handle.
- The effects of antibiotics other than gentamicin on the CanGaroo® Envelope have not been tested. The long-term effects of incorporating antibiotic treatments in the hydration of the CanGaroo® Envelope have not been assessed.

Caution: Federal Law (U.S.) restricts this device to sale by or on the order of a licensed medical practitioner.

POTENTIAL COMPLICATIONS

The following device-related complications are possible:

- Allergic reaction to ECM
 - Bleeding
 - Calcification
 - Fever
 - Fibrosis
 - Hematoma
 - Infection
 - Inflammation
 - Seroma
 - Undesired remodeling
- If the CanGaroo® Envelope is hydrated in gentamicin, complications associated with gentamicin are possible (neurotoxicity due to gentamicin, nephrotoxicity due to gentamicin, hypersensitivity to gentamicin)

STORAGE

This device must be stored in a clean, dry location at 10° C - 30° C.

STERILIZATION

This device has been sterilized with ethylene oxide gas.

REQUIRED MATERIALS

◆ Sterile dish (kidney dish or other bowl)

◆ Sterile forceps

◆ Suture

◆ Hydration fluid: a sufficient quantity of room temperature sterile water, sterile saline or sterile lactated Ringer's solution to completely immerse the CanGaroo® Envelope; at the surgeon's discretion, 20ml gentamicin 40mg/ml may be used to hydrate the Envelope¹

SUGGESTED INSTRUCTIONS FOR IMPLANTING THE CANGAROO® ENVELOPE

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. The CanGaroo® Envelope is MR safe in that it poses no known hazards in MR environments.

1. Select appropriate size CanGaroo® Envelope depending upon the size of the implantable electronic device that will be implanted. The CIED should fit snugly inside the CanGaroo® Envelope.
NOTE: The CanGaroo® Envelope cannot be used with implantable electronic devices that are larger than its internal dimensions.
NOTE: Ensure that the correct size Envelope is selected since sutured seam cannot be trimmed.
2. Inspect the CanGaroo® Envelope packaging/pouch for signs of damage. Do not use if package is opened or damaged.
3. Using aseptic technique, remove the inner Tyvek protective sleeve from the pouch, and place the inner Tyvek protective sleeve containing the CanGaroo® Envelope into the sterile field.
4. Open the Tyvek protective sleeve carefully and aseptically remove the CanGaroo® Envelope.
5. Hydrate the CanGaroo® Envelope by completely immersing it in a bowl of sterile water, sterile saline or sterile lactated Ringers solution for 1-2 minutes prior to use. If using 20ml gentamicin 40mg/ml for hydration, hydrate the CanGaroo® Envelope using aseptic technique by completely immersing it in the gentamicin for **2 minutes**.
6. Prepare the implantable electronic device as per manufacturer's instructions, making sure to secure the leads.
7. Slide the implantable electronic device into the pouch opening with lead wires emerging out of the envelope opening.
NOTE: Nonabsorbable or absorbable monofilament sutures can be used to tack the opening of the pouch to secure the implantable electronic device prior to implantation.
8. Place the implantable electronic device into the patient as per standard practice.


ADVERSE EVENT REPORTING

Any potential adverse incident involving the CanGaroo® Envelope should be reported immediately. Please report any device-related adverse events to Aziyo Biologics at +1 470-514-4080.

RETURN GOODS POLICY

For information on product returns and return authorization, contact Aziyo Biologics by calling +1 470-514-4080. All products returned to Aziyo Biologics must be accompanied by a Return Goods Authorization Number.

SYMBOLS AND THEIR EXPLANATIONS

	Use By
	Store at 10° C - 30° C
	Consult Instructions for Use
	Lot Number
	Catalog Number
	Sterilized using Ethylene Oxide
	Do Not Reuse
	Do Not Resterilize
	Do Not Use if the Package is Damaged or Open
	Non-pyrogenic
	Manufacturer
	Authorized Representative in the European Community
	MR Safe
	Contents

Aziyo® is a registered trademark of Aziyo Biologics, Inc.
© 2021 Aziyo Biologics, Inc. All rights reserved.

¹ Testing has shown that hydrating the CanGaroo® Envelope in gentamicin alone does not adversely affect the performance of the device. When hydrated in gentamicin prior to implantation, the CanGaroo® Envelope has been shown to reduce or inhibit bacterial colonization of the Envelope and the enclosed CIED in *in vitro* testing (*Methicillin-resistant Staphylococcus aureus* [MRSA], *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Pseudomonas aeruginosa*, and *Serratia marcescens*) and *in vivo* testing (*Methicillin-resistant Staphylococcus aureus* [MRSA], *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Pseudomonas aeruginosa*, and *Escherichia coli*) (Sohail MR et al. Preclinical evaluation of efficacy and pharmacokinetics of gentamicin containing extracellular-matrix envelope. Pacing Clin Electrophysiol. 2020;43(3):341-349). However, the claim of reduction of bacterial colonization of the Envelope and the enclosed CIED has not been established with human clinical data, nor has a clinical impact associated with this claim been demonstrated.

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Aziyo® är ett registrerat varumärke som tillhör Aziyo Biologics, Inc.
© 2021 Aziyo Biologics, Inc. Alla rättigheter förbehålles.